




## 2024 Enteral Feeding Devices Product Catalog





# TABLE OF CONTENTS

---

Section	Page #
ENFit Tip Syringes	6
Nasogastric (NG) Feeding Tubes	8
Pediatric (PU) & PVC NG Tubes	10
ENFit Pharmacy Items	12
Q Enteral Feeding Pump	14
Q Pump Feed Sets	18
Extension Feeding Sets	20
Gravity Feeding Sets	22
Gastrostomy Feeding Tubes (G Tubes)	24
Drainage Bags	26



# FEEDING LIFE

---

Innovative technology delivering  
a brighter future for tube feeders  
around the globe.




# Vesco

(vēscō)

Latin to English Translation  
to feed, to eat

## The Droplet

Symbolic of a drop of nutrition  
artistically flowing into carefully  
coiled tubing



## Vision

Worldwide leader delivering the  
broadest portfolio of the highest  
value enteral feeding devices



## Feeding Innovation

Our promise to deliver break-through  
enteral products to better serve our  
customers







## High-Quality Devices

---

- ISO 13485 Certified Quality Management System.
- CE Marked, FDA Registered & 510(k) Premarket Clearance.
- Made without natural rubber latex, DEHP, BPA or other phthalate plasticizers.



## Unbeatable Service

---

- Friendly, skilled and experienced customer care specialists.
- Full portfolio of ENFit products direct or through distributors.
- Responsive when and where you need us.



## Competitive Pricing

---

- Nimble, agile and driven to keep costs down and pass on savings.
- High-quality products without the premium prices.

# ENFit® Tip Syringes

## ENFit Tip Syringes

- Single Use Syringes compliant with ISO 80369-3 Standard
- MultiFil® Single Patient Use Syringes
- Sterile and Clean/Non-Sterile

## ENFit Standard Tip Syringes (5, 10, 20, 35 and 60mL)


- Single use enteral specific syringes for use in initial tube placement, flushing/hydration, bolus feeding administering medication and irrigation.
- Printed graduated length markings in mL.
- One cap included in each blister pack.
- Syringes are sold in clean/non-sterile and sterile configurations. Sterile devices are sterilized via ethyl-o-xylene (EO). Do not re-sterilize.

## ENFit Low Dose Tip (LDT) Syringes (0.5, 1 and 3mL)


- Designed to minimize process changes in pharmacy and bedside administration.
- Demonstrated dose accuracy administering smaller volumes compliant with GEDSA\*.
- Internal lumen consistent with traditional male syringe design reduces dead space and priming volume.
- One cap included in each blister pack.
- Device is sterilized (EO for ethyl-o-xylene). Do not re-sterilize.

## MultiFil ENFit Syringes (10 & 60mL)

- Single patient reusable syringe for up to 14 days if properly cleaned after each use. See instructions for additional information.
- Enteral specific syringes for use in initial tube placement, flushing/hydration, bolus feeding, medication administration and irrigation.
- Printed graduated length markings in mL.



ENFit  
Standard Tip




ENFit  
Low Dose Tip (LDT)

### Important Safety Information

- Use ENFit Tip Syringes only with ENFit feeding tube ports. Feeding tubes with different open female ports are not compatible with Vesco ENFit Tip Syringes.
- Made without natural rubber latex, Bisphenol A (BPA), DEHP or other phthalate plasticizers.
- Single-patient use.

ENFit is a federally registered trademark of GEDSA in multiple jurisdictions throughout the world. Used with permission. \* GEDSA August, 2016 Guidance Document

# ENFit Syringe Item Numbers



ENFit Tip Syringes						
Size (mL)	Single Use - Sterile		Clean/Non-Sterile		MultiFil	
	Item #	Units/Case	Item #	Units/Case	Item #	Units/Case
1 (LDT)	VED-601EO	800	N/A	N/A	N/A	N/A
3 (LDT)	VED-603EO	800	N/A	N/A	N/A	N/A
5	VED-605EO	800	VED-605	500	N/A	N/A
10	VED-610EO	800	VED-610	500	VED-610R	800
20	VED-620EO	600	VED-620	500	N/A	N/A
35	VED-635EO	400	VED-635	500	N/A	N/A
60	VED-660EO	400	VED-660	180	VED-660R	180
60 Thumb Ring	N/A	N/A	VED-661	180	N/A	N/A
60 Thumb Ring + 500mL Container	N/A	N/A	VED-665	180	N/A	N/A

ENFit Syringe Pharmacy Bulk Pack				
Size (mL)	Clear Item #	Units/ Case	Amber Item #	Units/ Case
0.5mL (LDT)	VED-600B	4000	VED-600AB	4000
1 (LDT)	VED-601B	4000	VED-601AB	4000
3 (LDT)	VED-603B	4000	VED-603AB	4000
5	VED-605B	2000	VED-605AB	1000
10	VED-610B	1000	VED-610AB	1000
20	VED-620B	1000	VED-620AB	1000
35	VED-635B	800	VED-635AB	800
60	VED-660B	400	VED-660AB	400


# Nasogastric Feeding Tubes

## Nasogastric Polyurethane (PU) Feeding Tubes

- Single or Dual ENFit Ports
- French Sizes: 5, 6, 8, 10, 12, 14 and 16
- Optional Stylet and/or Weighted Tip
- Lengths: 18, 36, 45 and 55 Inches

## Features

- Single or Dual ENFit Port for feeding and administering medication.
- Premium medical-grade kink-resistant polyurethane tubing.
- 4 side ports and distal end port improve flow and reduce likelihood of tube clogging.
- Water-activated lubricant aids in tube placement and allows for aseptic lubrication.
- Radiopaque for x-ray confirmation of tube placement.
- Centimeter graduated length markings.
- Available with and without flow through stylet for ease of placement and confirmation.
- Available non-weighted or non-magnetic tungsten weighted (closed tip).
- Not made with natural rubber latex, DEHP, Bisphenol A (BPA) or other phthalate plasticizers.



Single & Dual ENFit Port NG Tubes



## Important Safety Information

- Use only administration/feeding sets with ENFit connectors. Feeding sets with different connectors are not compatible with Vesco feeding tubes.
- Device is non-sterile.
- Tube has the potential to misconnect with small bore connectors of other healthcare applications.
- Single-patient use and Rx only.



# NG Tube Item Numbers

Single & Dual Port Polyurethane (PU) NG Tubes				
French Size	36" (91cm) Item #		45" (114cm) Item #	Units/ Case
	Single Port	Dual Port	Dual Port	
5	N/A	N/A	N/A	10
6	VED-19106	VED-29106	VED-21406	10
8	VED-19108	VED-29108	VED-21408	10
10	VED-19110	VED-29110	VED-21410	10
12	VED-19112	VED-29112	VED-21412	10
14	VED-19114	VED-29114	VED-21414	10
16	N/A	VED-29116	N/A	10

Single & Dual Port Polyurethane (PU) NG Tubes with STYLET				
French Size	36" (91cm) Item #		45" (114cm) Item #	Units/ Case
	Single Port	Dual Port	Dual Port	
6	N/A	VED-29106S	VED-21406S	10
8	VED-19108S	VED-29108S	VED-21408S	10
10	VED-19110S	VED-29110S	VED-21410S	10
12	VED-19112S	VED-29112S	VED-21412S	10
14	VED-19114S	VED-29114S	VED-21414S	10
16	N/A	VED-29116S	N/A	10

Dual Port Polyurethane (PU) NG Tubes, WEIGHTED with STYLET					
French Size	Weight (Grams)	36" (91cm) Item #	45" (114cm) Item #	55" (140cm) Item #	Units/ Case
6	0.55	VED-29106SW	VED-21406SW	N/A	10
8	0.75	VED-29108SW	VED-21408SW	VED-24008SW	10
10	1.33	VED-29110SW	VED-21410SW	VED-24010SW	10
12	2.07	VED-29112SW	VED-21412SW	VED-24012SW	10

# Pediatric (PU) & PVC NG Tubes

## Pediatric Polyurethane (PU) NG Tubes with ENFit®

- For use up to 30 days
- French Sizes: 4, 5, 6, 8, 10, and 12
- Lengths: 16, 24 and 35 Inches (40, 60 & 90cm)

### Features

- Clear and transparent tube for monitoring fluid content.
- Single ENFit port for feeding, flushing and administering medication.
- Indicated for long-term use up to 30 days.
- Premium medical-grade kink-resistant polyurethane tubing.
- 2 side ports improve flow and reduce likelihood of tube clogging.
- Radiopaque for x-ray confirmation of tube placement.
- Device is sterilized (EO for ethyl-o-xylene). Do not re-sterilize.
- Not made with natural rubber latex, DEHP, Bisphenol A (BPA) or other phthalate plasticizers.
- Centimeter graduated length markings.

## PVC NG Tubes with ENFit®

- French Sizes: 4, 5, 6, 8, 10, 12, 14, and 16
- Lengths: 16 and 35 Inches (40 & 90cm)

### Features

- Single ENFit port for feeding and administering medication.
- Indicated for short-term use, 24-48 hours.
- 2 side ports improve flow and reduce likelihood of tube clogging.
- Radiopaque for x-ray confirmation of tube placement.
- Not made with natural rubber latex, DEHP, Bisphenol A (BPA) or other phthalate plasticizers.
- Centimeter graduated length markings.



Pediatric & PVC  
NG Tubes



### Important Safety Information

- Use only administration/feeding sets with ENFit connectors. Feeding sets with different connectors are not compatible with Vesco feeding tubes.
- Devices are sterilized (EO for ethyl-o-xylene). Do not re-sterilize.
- Tube has the potential to disconnect with small bore connectors of other healthcare applications.
- Single-patient use.

# Pediatric & PVC NG Tube Item Numbers

Pediatric Polyurethane NG Tubes - Sterile				
French Size	16" (40cm) Item #	24" (60cm) Item #	35" (90cm) Item #	Units/Case
4	VED-84004EO	VED-86004EO	N/A	40
5	VED-84005EO	VED-86005EO	VED-89005EO	40
6	VED-84006EO	VED-86006EO	VED-89006EO	40
8	VED-84008EO	VED-86008EO	VED-89008EO	40
10	N/A	VED-86010EO	VED-89010EO	40
12	N/A	N/A	VED-89012EO	40

PVC NG TUBES - Sterile			
French Size	16" (40cm) Item #	35" (90cm) Item #	Units/Case
4	VED-94004EO	N/A	40
5	VED-94005EO	VED-99005EO	40
6	VED-94006EO	VED-99006EO	20
8	VED-94008EO	VED-99008EO	20
10	VED-94010EO	VED-99010EO	20
12	N/A	VED-99012EO	20
14	N/A	VED-99014EO	20
16	N/A	VED-99016EO	20



# ENFit Pharmacy Items

## ENFit Medicine Bottle Adapters (Multiple Sizes)

- Variety of sizes intended to meet most common medication bottles.
- Fits inside bottle opening to allow continued use of childproof caps.
- Ideal for aseptic technique.
- Made without natural rubber latex, DEHP, Bisphenol A (BPA) or other phthalate plasticizers.

## ENFit Transition Connectors

- Stepped Transition Connector for backwards compatibility to legacy feeding tubes.
- Available as a stand-alone item by the case.

## ENFit Tamper Evident Syringe Caps (Fits All Sizes)


- Provides security to syringe contents during transportation, prior to administration and enhances medication safety protocols.
- Made without natural rubber latex, DEHP, Bisphenol A (BPA) or other phthalate plasticizers.

## ENFit Syringe Tip Cap (Fits All Sizes)


- Pharmacy-friendly purple cap designed to aid in single handed placement.
- Ideal for aseptic technique and secure for transport.
- Made without natural rubber latex, DEHP, Bisphenol A (BPA) or other phthalate plasticizers.

## ENFit Medicine Straw (2, 4, and 6 inch)


- Connects to an ENFit Tip Syringe to easily facilitate the withdrawal of medication or feed.
- Ideal for aseptic technique.
- Device is sterilized (EO for ethyl-o-xylene). Do not re-sterilize.



Transition Connectors



ENFit Medicine  
Bottle Adapters



ENFit Medicine Straws

### Important Safety Information

- Use ENFit tip syringes only with ENFit feeding tube ports. Feeding tubes with different open female ports are not compatible with Vesco ENFit tip syringes.
- Single-patient use.




# Pharmacy Item Numbers



ENFit Medicine Bottle Adapters			
Size #	Inner Diameter (mm)	Item #	Units/Case
1	12 - 15	VED-690	100
2	14 - 16.5	VED-691	100
3	17 - 20.5	VED-692	100
4	19 - 23	VED-693	100
5	24 - 26	VED-694	100
6	27-30	VED-695	100
7	31-34	VED-696	100

Ancillary Items	Item #	Units/Case
Transition Connectors	VED-001	300
2 Inch (5cm) ENFit Medicine Straw - Sterile	VED-670EO	150
4 Inch (10cm) ENFit Medicine Straw - Sterile	VED-671EO	150
6 Inch (15cm) ENFit Medicine Straw - Sterile	VED-672EO	150
ENFit Syringe Screw-on-Cap	VED-697	4000
ENFit Syringe Pop-on-Cap (Dark)	VED-699P	4000
ENFit Syringe Pop-on-Cap (Lite)	VED-699PF	4000
ENFit Tamper Evident Cap	57-220	1000
ENFit Cap - Sterile	VED-699EO	4000



# Introducing You to The VESCO Q ENT

**New**

## SIMPLE

- + Easy Single-Handed Set Loading
- + Audible Click Ensures Correct Placement
- + Repeat Last Setting Enables Quick Programming & Consistent Feeds

## SMART

- + Interactive Touch Screen
- + Self-Guided Programming
- + Lock Screen Prevents Unwanted Interruptions



## VERSATILE

- ☒ Ambulatory for On-the-Go Use
- ☒ 1 to 600mL/hr up to 2000mL
- ☒ Bolus, Continuous & Intermittent Feeds
- ☒ Variable Alarm Volume & Brightness
- ☒ Auto Dimming Display



\*Do not use non-prescribed home-prepared, blenderized, or liquidized foods in the feed sets. Only use commercially available, prepared enteral feeding solutions and blenderized diets prescribed by the clinician.

# TERAL FEEDING PUMP



## STRONG

- + Feeds Thicker Formulas ( $\leq 2.5\text{Kcal}$ ) & Blended Diets\*
- + Drop, Scratch & Water Resistant
- + 24-Hour Battery Life



## RELIABLE

- ☒ +/- 5% Accuracy for All Feed Types
- ☒ Anti-Free Flow Set
- ☒ Easy to Clean Under Running Water (IP24)
- ☒ Precision Sensors for Optimal Audio & Visual Alarms

## QUIET

- ☒ No Loud Noises While Pumping
- ☒ Quieter than Competitive Pumps & Common Household Appliances



# Q Enteral Feeding Pump


**New**

## Vesco Q Enteral Feeding Pump

- Intermittent and Continuous Delivery of Enteral Fluids
- Intended for Adult & Pediatric Feeding Therapy in Hospital & Home Care Environments

### Q Pump Features:

- Includes region-specific power adapter & compact flexible pole clamp.
- Wide range of feed settings with power to feed thicker formulas  $\leq 2.5$  KCal & blended diets.
- User-friendly touchscreen interface with self-guided programming.
- Single-handed cassette loading & unloading.
- Ergonomic and portable with ambulatory functionality works in any orientation for on-the-go feeding.
- $\pm 5\%$  delivery accuracy of enteral fluids\*.
- Precision sensors for optimal audio and visual alarms.
- Water resistant and easy to clean with household or industrial cleaners.



Q Enteral Feeding Pump

### Important Safety Information

- Only use the pump to deliver enteral fluids at settings appropriate for the patient prescribed by a physician, registered dietitian, registered nurse, or other licensed practitioner.
- The battery capacity is an estimate that varies with use, flow rate, and run time. Charge before use.
- Do not clean pump while the AC Adapter/Charger is connected to the pump or while the pump is on. Do not connect AC Adapter/Charger to an outlet after cleaning until it is dry.
- Pump is not designed to deliver non-enteral fluids. Misuse may result in serious injury.
- The Vesco Q Pump does not require calibration. Pump should be returned to the recommended and approved service center for any service or repair requirements.



# Q Pump Item Numbers



## Product Specifications:

- Ambulatory:** Accuracy in upright orientation & on-the-go
- Weight:** ~ 1 lbs. (498 grams)
- Dimensions:** 3.75" x 6.5" x 2.0"  
(6.5 x 5.0 x 8.9cm)
- Display:** Capacitive Touch Screen
- Battery:** Lithium Ion  
24 Hour Life at 125mL/Hr  
6 Hour Charge Time
- Dose/Rate:** 1-600mL/Hr up to 2000mL
- Priming:** Manual or Auto
- Drop Resistance:** Any angle from 1 meter
- Classification:** Class II, B
- Rating:** IP24
- Warranty:** 24 Months



Vesco Q Pump	Item #	Units/Case
Q Pump - US with Type A Power Supply and Flexible Pole Clamp	VED-701	1
Q Pump - Type A Power Adapter - US 110	VED-731	1
Q Pump - Compact Flexible Pole Clamp	VED-738	1

# Q Pump Feed Sets

**New**

## Q Pump Feed Sets with ENFit

- ENPlus Spike & 40mm Screw Cap Closed System Sets
- 500mL & 1000mL PVC Bag Sets
- 1000mL Layered Laminate Bag Set


### Features:

#### ▪ Bag Sets

- Snap cap with rigid neck for easy filling and closure.
- Easy-to-read graduations and easy view translucent material.
- Tapered bag and exit port allows complete product drainage.

#### ▪ Closed System Sets

- Compliant to ISO 18250-3 Standard (ENPlus) for Enteral Reservoir Connectors.
- Screw cap allows for feeding through containers with a 40mm opening.
- Suspension bag for hanging containers on an IV Pole.
- Designed for single-handed cassette loading and unloading.
- Anti-Free Flow Silicone Tubing.
- Compliant to ISO 80369-3 Standard with ENFit Connectors.
- Not made with natural rubber latex, Bisphenol A, DEHP or other phthalate plasticizers.



### Important Safety Information

- Only use the pump to deliver enteral fluids at settings appropriate for the patient prescribed by a physician, registered dietitian, registered nurse, or other licensed practitioner.
- Only use Vesco Q Feed Sets with the Vesco Q Pump.
- Do not exceed the intended use time of the feed sets by programming a feeding regimen that exceeds 24 hours.
- Delivering non-enteral fluids may result in serious injury.

# Q Pump Set Item Numbers




1000mL PVC Bag Set



40mm  
Screw Cap Set

Vesco Q Pump & Sets	Item #	Units/Case
500mL Bag Q Pump Feed Set (PVC) with ENFit	VED-050	30
1000mL Bag Q Pump Feed Set (LL*) with ENFit	VED-051	30
40mm Screw Cap Q Pump Feed Set with ENFit	VED-052	30
ENPlus Spike Q Pump Feed Set with ENFit	VED-053	30
1000mL Bag Q Pump Feed Set (PVC) with ENFit	VED-054	30

\*LL = Layered Laminate




# Extension Feeding Sets

## Universal Low-Profile Feeding Tube Extension Sets

- Bolus Feed Single ENFit Port with Straight or Right Angle available in 12" (30cm) and 24" (60cm)
- Continuous Feed Dual ENFit Port available in 12" (30cm) & 24" (60cm)

### Features

- Compatible with MIC-KEY® and Mini ONE® and all other low-profile gastrostomy feeding tubes.
- Dual ENFit Ports with 90 degree connector allows for continuous feeding, easy flushing and administration of liquid medication without disconnecting set.
- Single ENFit port & straight connector for bolus feeding, flushing and medication delivery.
- Not made with natural rubber latex, Bisphenol A (BPA) or other phthalate plasticizers.




Low Profile  
Extension Sets

## ENFit to ENFit Extension Sets

- Single & Bifurcated Feeding Ports
- Lengths: 20" (50cm), 39" (100cm) & 59" (150cm)

### Features

- ENFit male Port(s) to ENFit female provide longer length tubing to connect to ENFit compatible feeding devices.
- Available with in-line ENFit Y-Port for continuous feeding, easy flushing and administration of liquid medication without disconnecting set.
- Made without natural rubber latex, Bisphenol A, DEHP or other phthalate plasticizers.



ENFit to ENFit  
Extension Sets


### Important Safety Information

- Use only administration/feeding sets with ENFit connectors. Feeding sets with different connectors are not compatible with Vesco feeding tubes.
- See individual product labeling for directions and important safety information.
- Single-patient use.

The Avanos MIC-KEY and AMT Mini ONE are not registered trademarks of Vesco Medical.  
All Rights Reserved.



# Extension Set Item Numbers




## Universal Low-Profile Extension Sets

Set Description	12" (30cm) Item #	24" (60cm) Item #	Units/Case
Continuous - Dual Port, Right Angle	VED-1000	VED-1001	10
Bolus - Single Port, Straight	VED-1002	VED-1003	10
Bolus - Single Port, Right Angle	VED-1004	VED-1005	10

**New**

## ENFit to ENFit Extension Sets - Sterile

Set Description	Item #	Units/Case
20" (50cm) Bifurcated ENFit Ports to ENFit Extension Set	VED-1020EO	10
39" (100cm) Single ENFit Port to ENFit Extension Set	VED-1040EO	10
59" (150cm) Co-Extruded Single ENFit Port to ENFit Extension Set	VED-1059EO	10
59" (150cm) Bifurcated ENFit Ports to ENFit Extension Set	VED-1060EO	10




# Gravity Feeding Sets

## Gravity Feeding Sets with ENFit

- 500 and 1000mL Gravity Bag Sets
- ENPlus Spike and 40mm Screw Cap Closed System Gravity Sets
- Large Bore 1000mL Gravity Bag Set

## Features

- Gravity Bag Sets
  - Snap cap with rigid neck for easy filling and closure.
  - Easy-to-read graduations and easy-to-view translucent material.
- **Closed System Sets**
  - Screw cap allows for flexibility in feeding containers with a 40mm opening.
  - Convenient suspension bag allows for use with smaller containers (baby bottles) for easy hanging.
- **Large Bore Gravity Sets**
  - Designed to feed thicker formulas faster with 30° tapered bag for quicker, more complete feed.
  - Snap cap with rigid neck for easy filling and closure.
- Roller clamp and 14 drip chamber for easy priming and consistent flow.
- Not made with natural rubber latex, DEHP, Bisphenol A (BPA) or other phthalate plasticizers.
- Device intended to deliver liquid nutrition formulas or water to a patient's feeding tube.



1000mL  
Gravity Feed Set

## Important Safety Information


- To avoid risk of a tubing misconnection, ONLY connect the distal connector to an ENTERAL feeding tube. This device is compatible with an ENFit female connector on the patient's enteral feeding tube.
- Not for IV use.
- This device should be replaced at least every 24 hours.
- Do not store for prolonged periods at extreme temperatures.
- See individual product labeling for directions and important safety information.
- Single-patient use.

# Gravity Set Item Numbers



Gravity Feed Sets	Item #	Units/Case
40mm Screw Cap Set with ENFit	VED-022	30
ENPlus Spike with ENFit	VED-026	30
500mL Bag Set (LL)* with ENFit	VED-040	30
500mL TOP FILL Bag Set (PVC) with ENFit	VED-042	30
1000mL Bag Set (LL)* with ENFit	VED-034	30
1000mL TOP FILL Bag Set (PVC) with ENFit	VED-046	30
1000mL LARGE BORE Bag Set with ENFit	VED-049	30

\*LL = Layered Laminate



# Gastrostomy Feeding Tubes

## Dual ENFit Port Replacement Gastrostomy Feeding Tubes

- French Sizes: 14, 16, 18, 20, 22 and 24

### Features

- Dual ENFit Port allows for easy flushing and administration of liquid medication without disconnecting set.
- Medical-grade silicone tubing with clear shaft and graduation marks in 1cm increments.
- Rounded skin disc for retention.
- Gastric balloon for internal retention with rounded tip, open distal end and radiopaque tip marker.
- Device is sterilized (EO for ethyl-o-xylene). Do not re-sterilize.
- Not made with natural rubber latex, Bisphenol A (BPA) or other phthalate plasticizers.



Dual Port  
Gastrostomy Tube



### Important Safety Information

- Use only administration/feeding sets with ENFit connectors.
- Single-patient use.

# G-Tube Item Numbers



Dual ENFit Port Gastrostomy Tubes - Sterile			
French Size	Balloon Size (mL)	Item #	Units/Case
14	5	VED-214	5
16	5	VED-216	5
18	20	VED-218	5
20	20	VED-220	5
22	20	VED-222	5
24	20	VED-224	5




# Drainage Bags

## Enteral Drainage Bags with ENFit

- 250, 500, 1000 and 1500mL Drainage Bag Sets
- 500, 1000, and 1500mL Drainage Bag with 150cm proximal tubing

## Features

- Designed specifically for drainage in enteral feeding systems with male and female ENFit connectors compliant with ISO 80369-3 Small Bore Connector Standards.
- Large, clear tubing allows for visualization of gastric contents and facilitates gastric decompression or drainage.
- Eliminates the need for non-enteral drainage products such as urinary bags.
- Tethered cap & thumb clamp.
- Larger sizes (500, 1000, & 1500mL) equipped with 60 inches (150 centimeters) of extension tubing.
- Smaller sizes (250mL) have 30 inches (75 centimeters) of proximal tubing.




## Important Safety Information

- Use only with enteral feeding tubes with ENFit connectors.
- Device is sterilized (EO for Ethyl-o-xylene). Do not re-sterilize.
- Single-patient use.

# Drainage Bag Item Numbers



Drainage Bag Sets - Sterile	Item #	Units/Case
250mL PVC ENFit Drainage Bag with 30" (75cm) tubing	VED-009EO	30
500mL PVC ENFit Drainage Bag with 60" (150cm) tubing	VED-010EO	30
1000mL PVC ENFit Drainage Bag with 60" (150cm) tubing	VED-012EO	30
1500mL PVC ENFit Drainage Bag with 60" (150cm) tubing	VED-015EO	30



# Vesco Medical

Tube feeding experts delivering high quality products with unbeatable service at a competitive price. Our portfolio is All ENFit All The Time!

## Group Purchasing Organization Agreements:

HealthTrust (HT31810) MHA (MS60238N)

Premier (PP-DI-1965) Vizient (MS6216)

[www.VescoMedical.com](http://www.VescoMedical.com)

[info@vescomedical.com](mailto:info@vescomedical.com)

(855) 239-0131

**GEDSA**

ENFit is a registered trademark of GEDSA.  
©2024 Vesco Medical, LLC MKT 1128-8-05

 stay connected®

For more information about ENFit,  
visit [www.stayconnected.org](http://www.stayconnected.org)

 **VESCO**  
MEDICAL  
**Feeding Innovation**  
60 Collegeview Rd.  
Westerville, OH 43081