

# AMSafe® IV Pump Administration Sets






For use with Baxter Spectrum®, Flo-Gard® 6201 and 6301 series pumps.

## Features



- **Sure-Lok®** Needle-Free connector and Pre-Pierced Y-injection site options
- 0.22 micron filter option available
- TUBETECH® technology provides +/- 5% accurate fluid delivery over 24 hour timeframe
- 510K approved IV pump sets for use with Baxter Spectrum and Flo-Gard® 6201 and 6301 series pumps
- DEHP-Free, Latex-Free, and Non-Pyrogenic




## AMSafe® IV Pump Administration Sets with **Sure-Lok®** Needle-Free Connector

Item Code	Drops/ mL	Length	Description	Approx Priming Volume	Qty/Unit
<b>FP092</b>	10	92 in. (234 cm)	Roller Clamp, Slide Clamp, Rotating Male Luer Lock with Retractable Collar	14 mL	50/Case
					
<b>FP092Y1</b>	10	92 in. (234 cm)	1 Needle-Free Valve, Roller Clamp, 2 Slide Clamps, Rotating Male Luer Lock with Retractable Collar	14 mL	50/Case
					
<b>FP100BY2</b>	10	100 in. (254 cm)	Back check Valve, 2 Needle-Free Valves, Roller Clamp, 2 Slide Clamps, Rotating Male Luer Lock with Retractable Collar	16 mL	50/Case
					
<b>FP109BY3</b>	10	109 in. (277 cm)	Back check Valve, 3 Needle-Free Valves, Roller Clamp, 2 Slide Clamps, Rotating Male Luer Lock with Retractable Collar	18 mL	50/Case
					
<b>FP169BY3</b>	60	109 in. (277 cm)	Back check Valve, 3 Needle-Free Valves, Roller Clamp, 2 Slide Clamps, Rotating Male Luer Lock with Retractable Collar	18 mL	50/Case
					

## AMSafe® IV Pump Administration Sets with **Sure-Lok®** Needle-Free Connector and Filter


Item Code	Drops/ mL	Length	Description	Approx Priming Volume	Qty/Unit
<b>FP106FY2</b>	10	106 in. (269 cm)	Back check Valve, 2 Needle-Free Valves, Roller Clamp, 2 Slide Clamps, 0.22 Micron Filter, Rotating Male Luer Lock with Retractable Collar	22 mL	50/Case
					
<b>FP100FY1</b>	10	100 in. (254 cm)	1 Needle-Free Valve, Roller Clamp, 2 Slide Clamps, 0.22 Micron Filter, Rotating Male Luer Lock with Retractable Collar	22 mL	50/Case
					

For use with Baxter Spectrum®, Flo-Gard® 6201 and 6301 series pumps.




### IV Pump Administration Sets with Pre-Pierced Y-Injection Site

Item Code	Drops/ mL	Length	Description	Approx Priming Volume	Qty/Unit
<b>FP100BI2</b>	10	100 in. (254 cm)	Back check Valve, 2 Pre-Pierced Injection Y-Sites, Roller Clamp, 2 Slide Clamps, Rotating Male Luer Lock with Retractable Collar	16 mL	50/Case




<b>FP109BI3</b>	10	109 in. (277 cm)	Back check Valve, 3 Pre-Pierced Injection Y-Sites, Roller Clamp, 2 Slide Clamps, Rotating Male Luer Lock with Retractable Collar	18 mL	50/Case
-----------------	----	---------------------	--	-------	---------



### IV Pump Administration Sets with Pre-Pierced Y-Injection Site and Filter

Item Code	Drops/ mL	Length	Description	Approx Priming Volume	Qty/Unit
<b>FP106FI2</b>	10	106 in. (269 cm)	Back check Valve, 2 Pre-Pierced Injection Y-Sites, Roller Clamp, 2 Slide Clamps, 0.22 Micron Filter, Rotating Male Luer lock with Retractable Collar	22 mL	50/Case



To view our complete product line, visit us at [www.amsino.com](http://www.amsino.com)



#### Amsino International, Inc.

708 Corporate Center Drive, Pomona, CA 91768  
800.632.6746 • fax: 909.626.3888 • [info@amsino.com](mailto:info@amsino.com)

Product(s) described may not be licensed or available for sale in all countries. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. AMSafe and Sure-Lok are registered trademarks of Amsino International, Inc. © Copyright 2019, Amsino International, Inc. All rights reserved.